

Product datasheet for **RR208445L4V**

Syng2 (NM_053553) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Syng2 (NM_053553) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Syng2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053553
ORF Size:	702 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208445).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053553.1 , NP_446005.1
RefSeq Size:	1108 bp
RefSeq ORF:	705 bp
Locus ID:	89815
UniProt ID:	O54980
Cytogenetics:	10q32.2
Gene Summary:	ubiquitously expressed synaptogyrin isoform; undergoes tyrosine phosphorylation by pp60c-src; may play a role in membrane trafficking [RGD, Feb 2006]



[View online »](#)