

## Product datasheet for **RR208448L4V**

### Srsf6 (NM\_001014185) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Srsf6 (NM_001014185) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Srsf6
Synonyms:	Sfrs6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001014185
ORF Size:	1017 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208448).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001014185.2</a> , <a href="#">NP_001014207.2</a>
RefSeq Size:	3513 bp
RefSeq ORF:	1020 bp
Locus ID:	362264
UniProt ID:	<a href="#">G3V6S8</a>
Cytogenetics:	3q42



[View online »](#)

**Gene Summary:**

Plays a role in constitutive splicing and modulates the selection of alternative splice sites. Plays a role in the alternative splicing of MAPT/Tau exon 10. Binds to alternative exons of TNC pre-mRNA and promotes the expression of alternatively spliced TNC. Plays a role in wound healing and in the regulation of keratinocyte differentiation and proliferation via its role in alternative splicing (By similarity).[UniProtKB/Swiss-Prot Function]