

Product datasheet for **RR208496L4V**

Adra1d (NM_024483) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adra1d (NM_024483) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adra1d
Synonyms:	Adrd1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024483
ORF Size:	1683 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208496).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024483.1 , NP_077809.1
RefSeq Size:	2939 bp
RefSeq ORF:	1686 bp
Locus ID:	29413
UniProt ID:	P23944
Cytogenetics:	3q36
Gene Summary:	adrenergic receptor showing high affinity for the antagonist WB4101 and agonists phenylephrine and methoxamine [RGD, Feb 2006]



[View online »](#)