

Product datasheet for **RR208520**

Cbx8 (NM_001034078) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cbx8 (NM_001034078) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cbx8
Synonyms: MGC116381
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR208520 representing NM_001034078
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGCTCTCGGCGTGGGGAGCGGGTGTTCGCGGCCGAAGCCCTCTGAAGCGGCGCATTTCGAAAG
GACGCATGGAATATCTCGTGAAATGGAAGGGCTGGTCGAGAAGTACAGCACGTGGGAGCCTGAAGAAA
TATTCTTGATGCTCGCCTCCTTGCAGCCTTTGAAGAAAGGAACGGGAGATGGAGCTCTATGGTCCCAA
AAGCGAGGACCCAAGCCTAAAACCTTCTCTCAAGGCCAGGCCAAGGCAAGGCCAAAACCTATGAAT
TCAGAAGCGACTCTACCAGAGGCATCCGGATCCCTTACCCAGGCCGCTCTCCCCAGGATTTGGCATCTAC
TTCTAGGGCCCGAGAGGGCCTTCGGAACACGGGGCTACCCCCACCGGGAGCAGTACCAGTACCTGCAGG
GCAGACCCACCTCGGGACAGAGACCGAGACCGGGATCGAGAACGGGACAGGGGTACCAGCCGTGTAGATG
AGAAGCCCAGCTCACCGGGGACAGCTCCAAGAAACGAGGACCCAAGCCCAGGAAGGAGCTCCTAGACCC
TTCACAGAGACCTTTGGGAGAGCCAGTGATGGCCTTGGAGAATACCTCAAAGGCAGGAAGCTGGATGAG
ACCCCTCTGGGATAGGAAAGTTCTCAGCTGGCCACAGCGTGATCCAGCTTGTCTGGAGGCAGGACTCAG
ACCTGGTACAATACGGTGTGACCAGTCTAGCTCGGCGGAGGCCCTCGGCAAGTTGGCTCGGACACCTT
TCCAGCCAGGGTGATAAAGCACCGGGCTGCTTTACTGGAGGCCAAGGTCAAGGTGCCCTGGACCCCTGGT
GGGACCAGGGTCCGACATAGTTCAGGCACCCCGGGCTCAGTGGGAAGCCTGTATCGGGACATGGGGCTC
AAGGGGAAGGCCCTCCCTCATCGCCAGGATCCCAGTGGCCAGAATCCTGGGGGACCCGGAGGAAGATC
CTGGAGCCCTCTCTGACTAACCTGGAGAAGGTGGTCGTACAGACGTGACCTCAAACCTTTTGGACCGTC
ACCATTAAGGAGAGTAACACGGACCAAGGCTTCTTTAAGGAGAAAAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR208520 representing NM_001034078
 Red=Cloning site Green=Tags(s)


MELSAVGERVFAAEALLKRRIRKGRMEYLVKWKGWSQKYSTWEPEENILDARLLAAFEEREREMELYGPK
 KRGPKPKTFLLKAQAKAKAKTYEFRSDSTRGIRIPYGRSPQDLASTSRAREGLRNTGLPPPGSSTSTCR
 ADPPRRDRDRDRERDRGTSRVDEKPPSSPGDSSKKRGPKPRKELLDPSSQRPLGEPDGLGEYLKGRKLDE
 TPSGIGKF SAGHSVIQLARRQSDLVQYGVTSPPSSAEASGKLAADTFPARVIKHRAALLEAKGQGALDPG
 GTRVRHSSGTPGSVGSLYRDMGAQGRPSLIARIPVARILGDPEEESWSPSLTNLEKVVVTDVTSNFLT
 VTIKESNTDQGFKEKR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001034078

ORF Size:	1098 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001034078.1 , NP_001029250.1
RefSeq Size:	1383 bp
RefSeq ORF:	1101 bp
Locus ID:	303731
Cytogenetics:	10q32.3
MW:	40.4 kDa