

Product datasheet for **RR208582L4V**

Lat2 (NM_173840) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lat2 (NM_173840) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lat2
Synonyms:	Ntal; Wbscr5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_173840
ORF Size:	612 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208582).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_173840.1 , NP_776212.1
RefSeq Size:	615 bp
RefSeq ORF:	615 bp
Locus ID:	317676
UniProt ID:	Q8CGL2
Cytogenetics:	12q12
Gene Summary:	human homolog maps to the microdeletion region in Williams syndrome and may act as an adaptor for non T-cell immune response [RGD, Feb 2006]



[View online »](#)