

Product datasheet for **RR208591**

Ppil2 (NM_001017383) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ppil2 (NM_001017383) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ppil2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR208591 representing NM_001017383
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGGAAGCGACAGCACCAGAAGGACAAGATGTACATCACCTGTGCAGAGTACACTCATTCTATGGTG
 GCAGGAAGCCAGATATCACACAGACAAGTTTTCGCCGCTTACCTTTTGACCATTGCAGTCTCTCTCCA
 GCCTTTTGTCTACCCAGTCTGCACCCCGGAAGGTGTCGTCTTCGACTTGCTGAACATTGTCCCATGGCTT
 AAGAAGTATGGGACCGATCCGAGCACTGGAGAGAACTTGATGGGAAGTCCCTGATTAAGCTGAACCTTTG
 CAAAGAACAGTGAAGGGCAGTACCCTGTCCAGTGTACTCCGTGTTCACTGACAACCCACATCGT
 GGCCATCAGGACAACGGCAATGTCTACACCTATGAGGCAGTGGAGCAGCTAAACATCAAGGCCAAGAAC
 TTGAGGGACCTGTTAACTGATGAGCCCTTTCCAGGCAAGACATCATACCCTGCAGGACCCCAACATC
 TGGACAAATTTAACGTCAGCAGCTTCTCCATGTAAGAATAACATGCGAATGATAGATCCAGATGAGGA
 AAAGGCCAAACAAGCCATCTTATTATTTGAAAATACAAATGCGGAGACCCGAGAGACGCTACAGGAG
 CTCTACAAAGAGTTCAAAGGAGATGAGATTTTAGCAGCTACCATGAAGCCACCTGAGAAGAAGAAGGTGG
 ACCAACTGAATGCCGCCACTATTCCACAGGGAAGGTCAGCGCGTCTTACCTCTACCGCCATGGTGCC
 TGAGACCACGCATGAAGCAGCTGTATTGATGAAGATGCTCTGCGCTACCAGTTTGTGAAGAAGAAGGGC
 TATGTGAGGCTACACACCAACAAGGGTGACCTCAACCTCGAGCTGCACTGTGACCTGACACCAAAAACCT
 GTGAAAACCTCATCAAGCTCTGCAAGAAACAGTATTATGATGGTACCATCTTTCACAGATCCATCAGGAA
 CTTTGTGATCCAGGGTGGTGACCCTACTGGTACAGGCACAGGAGGAGTCACTGAGGCAAGCCTTTC
 AAAGACGAGTCCGTCCTCAATCTTTCACACACAGGCCGTTGGGTGCTCAGCATGGCCAATTCAGGGCCCA
 ACACCAACAAGTCTCAGTCTTTCATCACCTTCCGCTGTGCTTACCTGGACAAGAAGCACACCATCTT
 TGGACGGGTTGTCGGGGGCTTTGACACGCTGACAGCTATGGAGAATGTGGAGAGTGACCCCAAACTGAC
 CGCCCTAAGGAAGAAGTGCATCTGTACCACCACAGTGTGTTGTGGATCCCTATGAGGAGCCGACGCC
 AGATTGCCAGGAACGGAAGAAGACACAGCACCAAGTGGATCCAGAGGCCAAGGCCAAGCAAGCAAGCCAGCC
 CCAGCCTGGAACCCAGGGCCCCAGACATACCGCCAGGGGGTGGGCAAGTACATCAACCCTGCAGCCACG
 AAACGAGCAGCAGAGGAAGAACCATCCACCAGCACAGCTATCCCCACAGCCAAGAAAAGGCCAGTCGGG
 GCTTTGGGGACTTCAGCTCCTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR208591 representing NM_001017383
 Red=Cloning site Green=Tags(s)


MGKRQHQQDKMYITCAEYTHFYGGRKPDITQTSFRRLPFDHCSLSLQPFVYVCTPEGVVFLLNIVPWL
 KKYGTDPSSTGEKLDGKSLIKLNFKNSEGQYHCPVLYSVFTDNTHIVAIRTGNVYTYEAVEQLNIKAKN
 LRDLLTDEPFSRQDIITLQDPTNLDFNVSSFFHVKNMRMIDPDEEKAKQDPSYLLKNTNAETRETLQE
 LYKEFKGDEILAAATMKPPEKKKVDQLNAAHYSTGKVSASFSTAMVPETTHEAAVIDEDVLRYPVKKKG
 YVRLHTNKGDLNLELHCDLTPKTCENFIKLCCKQYYDGTIFHRSIRNFVIQGGDPTGTGTGGSYWGKPF
 KDEFRPNLSHTGRGVLSMANSGPNTNKSQFFITFRSCAYLDDKHTIFGRVVGGFDLTAENVESDPKTD
 RPKEEVRICTTTVFVDPYEEADAQIAQERKKTQHVDPEAKAKTNQPPGNQGPQYRQGVGKYINPAAT
 KRAAEEEPSTSTAIPTAKKRPSRFGDFSSW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001017383

ORF Size: 1563 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001017383.1](#), [NP_001017383.1](#)

RefSeq Size: 1801 bp


RefSeq ORF: 1566 bp

Locus ID: 360746

Cytogenetics: 11q23

MW: 59 kDa

Product images:



Circular map for RR208591