

Product datasheet for **RR208631**

S1pr4 (NM_001108075) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	S1pr4 (NM_001108075) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	S1pr4
Synonyms:	Edg6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR208631 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCAAACAACAGGAAACCCAGGCCAGGCTGTGACTCCACCCTGCCGGGGCCGCTCTCAGTCTATCCC
CGGGGAAGGCCATGAACATCAGTACCTGGTCCACGCTAGTGACCCAGAGTCTGCCACCGACTGGCGGC
TAGCGGTACAGCCTGCTCATCGTCTGCACTACAATCACAGCGGCAGGCTGGCTAGCCGTGGGGCCCT
GAGGACAAATGGCGGGCTGGGATGCTGCGGGGGCCGTCGGTGGCAGCGGGCTGCCTGGTGGTCTGGAAA
ACGCCATGGTGTGGCTGCCATCGCCATCCACATGCGCTCACGGCGCTGGGTGTACTACTGCCTCCTGAA
CATCACCCCTGAGCGACCTGCTCACAGGCCCTGGCCTACGTGGTCAACGTGCTGCTGTAGGAACTCGCACC
TTCCAGCTGTACCTGTGCACTGGTTCCTGCGGGAGGGCCTGCTCTTCATGGCCCTGGCTGCGTCCACCT
TCAGTCTCCTCTTCACAGCCGGGAGCGCTTCGCCACCATGGTGCGGGTGGCTGAGAGCGGGGCCACCAA
GACCAGCCGTGTGTATGGCTGCATCGGTCTGTGCTGGCTGCTGGCAGCTACCCCTGGGCTGCTGCCCTG
CTGGGCTGGAAGTGTGTGCGCCTCCAGCGCTGCTCTAGCCTGCTGCCCTCTACTCCAAGGGCTATG
TGCTCTTTGTGTGGTGGTCTTCGCCCTAATCCTAGTACTATCCTGAGCCTCTACGGGGCCATCTTTAG
GGTGGTCCGAGCCAACGGGCAGAAGTCCCCGCTCCTCCTGCCCGCCAAGTCCCGCAGGCTACTCAAC
ACCGTCTGATGATCTTGGTGGCTTTTGTGGTGTGCTGGGGTCCCCTGTTTGGCCTGCTCCTGGCCGACA
TCTTTGGATCTAATGCTGGGCCAGGAGTACCTGCGCGCATGGACTGGATCCTGGCCCTAGCTGTGCT
CAACTCAGCCATCAATCCTCTCATCTATTCTCCGAGCCGTGAGGTGCAGCACGCTGTGCTGACCTTC
CTGTGCTGCGGCTGCCTCAGGTTAGGCCTGAGAGGCCCTGGAGACTGCCTGACCCGGATCACCGAGGCC
ACTCTGGGGCATCCACCACTGACAGCTCGCTGAGGCCAGGGACAGTTTTTCGGACTTCGAGGTCACTCAG
CTTCAAGATGCGAGAGCCGCTGTCCAGCGTTTCCAGCATCCGAGCGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR208631 protein sequence
 Red=Cloning site Green=Tags(s)

MSNNRKRPRPGCDSTLPGAALSLSPGKAMNISTWSTLVTPESCHRLAASGHSLILVLYHNSGRLASRGGP
 EDNGGLGMLRGPVSAAGCLVLENAMVLAATAIHMRSRRWVYVYCLLNITLSDLLTGLAYVVNVLLSGTRT
 FQLSPVHWFLREGLLFMALAASTFSLFTAGERFATMVRVAESGATKTSRVYGCIGLCWLLAATLGLLPL
 LGWNCVCAFQRCSLLPLYSKGYVLFVVFVALILVTILSLYGAIFRVVRANGQKSPRPPARRKSRLLN
 TVLMILVAFVVCWGLFGLLLADIFGSNVWAQEYLRGMDWILALAVLNSAINPLIYSFRSREVQHAVLTF
 LCCGCLRLGLRPGDCLTRITEAHSGASTDSSLRPRDSFRTSRSLFSKMRPLSSVSSIRSA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001108075

ORF Size: 1239 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001108075.1, NP_001101545.1

RefSeq Size: 2277 bp


RefSeq ORF: 1242 bp

Locus ID: 314649

Cytogenetics: 7q11

MW: 44.9 kDa

Product images:



Circular map for RR208631