

Product datasheet for **RR208644L3V**

Plpp3 (NM_138905) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Plpp3 (NM_138905) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Plpp3
Synonyms:	Dri42; Ppap2b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138905
ORF Size:	936 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208644).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138905.2 , NP_620260.2
RefSeq Size:	3105 bp
RefSeq ORF:	939 bp
Locus ID:	192270
Cytogenetics:	5q34
Gene Summary:	an endoplasmic reticulum resident protein; upregulated during epithelial differentiation [RGD, Feb 2006]



[View online »](#)