

Product datasheet for **RR208663L3V**

Eloa (NM_017103) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eloa (NM_017103) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Eloa
Synonyms:	El; Tceb3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017103
ORF Size:	2319 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208663).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017103.2 , NP_058799.2
RefSeq Size:	2745 bp
RefSeq ORF:	2322 bp
Locus ID:	25562
Cytogenetics:	5q36
Gene Summary:	subunit A of the transcription factor B (SIII) complex (also known as elongin); frees RNA Polymerase II from transient pausing at template-encoded arresting sites [RGD, Feb 2006]



[View online »](#)