

Product datasheet for **RR208667**

Gabrr2 (NM_017292) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gabrr2 (NM_017292) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gabrr2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR208667 representing NM_017292 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGC**C

ATGCCTTATTTATGAGACTCGCTTTGTTCTTGTGTTTGCCTGATGGCTCTTGTGGAGAGCAGAAAACCA
GGAGGAAGAGATGGACAGGGCACCTGGAACATCCAAGCCAAGTCACTTGTATAAAAAGAACCTCGATGT
GACCAAGATCCGGACAGGAAAACCTCGGCCTGTCTCAGGGTAGAGGACCATGACTTCACCATGAGGCCA
GCCTTTGGAGGCCCGCCATCCCTGTGGGTGTGGACGTGCAGGTGGAGAGTCTGGACAGCATCTCCGAGG
TGGACATGGACTTCACCATGACCCTGTACCTGAGGCATTACTGGAGGGATGAGAGGCTGGCCTTTCCAG
CAGCAGCAACAGGAGCATGACCTTCGATGGCAGGCTGGTGAAGAAGATCTGGGTCCTGACGTCTTCTTT
GTTCACTCCAAAAGGTCAATTTACCCACGACACCACCACCGACAACATCATGCTGCGAGTGTTCGCGATG
GCCACGTGCTGTACAGCATGAGGATTACCGTCACTGCCATGTGCAACATGGACTTCAGCCACTTCCCCT
GGACTCCCAGACCTGCTCTTTGGAGCTAGAGAGCTACGCGTACACGGACGAGGATCTGATGCTGTACTGG
AAGAACGGGGATGAGTCCCTGAAGACGGACGAGAAGATCTCCTTGTCCCAGTTTCTGATCCAGAAGTTCC
ACAGACTTCCAGACTGGCCTTCTACAGCAGCACGGGCTGGTACAACCGTCTGTACATCAACTTCACTCT
GCGACGTACATCTTCTTCTTCTGCTCCAAACCTACTTCTGCCACCCTCATGGTCATGCTGTCTCTGG
GTGTCTTCTGGATCGACCACAGAGCTGTGCCTGCTAGAGTTTCACTGGGCATCATGACGGTGTGACCA
TGTCCACCATCATCACGGGCGTGAACGCTCCATGCCTCGAGTGTCTACATCCGGGCCGTGGACATCTA
CCTCTGGGTGACGCTTTGTGTTCTGTTTCTCTCAGTGTGGAGTACGCAGCTGTCAACTACCTGACCACA
GTGCAGGAACAGAAGGAACGGAAGCTTCGGGACAAATTTCCATGCACGTGTGGAATGCTTCATTCAAGAA
CCATGACGCTGGATGGAAGCTACAGTGAATCTGAGGCCAACAGCCTGGCTGGCTACCCAAGAAGCCACAT
TCTTCCAGAAGAAGAAAGGCAAGACAAAATAGTGGTTTCAATTTGGCCCTCAACAGTGAATTGACTTCATCC
AGAAAGAAAGGCTTCTCAAAGGCCAAATGGGACTCTATATCTCCAGAACACCCATGCCATTGACAAAT
ACTCCAGTTGATATTTCTGCCTTCTACATAGTTTTCAACCTAATTTATTGGTCAGTGTCTTCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR208667 representing NM_017292
 Red=Cloning site Green=Tags(s)

MPYFMRLALFLFCLMALVESRKPRRKRWTGHLETSKPSHLYKKNLDVTKIRTGKPRLLRVEDHDFTMPR
 AFGGPAIPVGVQVESLDSISEVDMDFMTLYLRHYWRDERLAFSSSNRSMTFDGRLVKKIWVPDVFF
 VHSKRSTHDTTDDNIMLRVFPDGHVLYSMRITVTAMCNMDFSHFPLDSQTSLELESYAYTDEDLMLYW
 KNGDESLKTDEKISLSQFLIQKFHTTSRLAFYSSTGWYNRLYINFTLRRHIFFFLLQTYFPATLMVMLS
 VSFWIDHRAVPARVSLGIMTVLTMSTIITGVNASMPRVSYIRAVDIYLWVSFVFLSVLEYAAVNYLTT
 VQEQKERKLRDKFPCTCGMLHSRTMTLDGSYSESEANSLAGYPRSHILPEEERQDKIVVHLALNSELTS
 RKKGLLKGMGLYIFQNTHAIDKYSRLIFPAFYIVFNLIYWSVFS


SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN:	NM_017292
ORF Size:	1395 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017292.2 , NP_058988.1
RefSeq Size:	1757 bp
RefSeq ORF:	1398 bp
Locus ID:	29695
UniProt ID:	P47742
Cytogenetics:	5q21
MW:	54.3 kDa
Gene Summary:	may form bicuculline-insensitive GABA channels in the retina [RGD, Feb 2006]