

Product datasheet for **RR208668L4V**

Pla1a (NM_138882) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pla1a (NM_138882) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pla1a
Synonyms:	Ps-pla1; Pspla1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138882
ORF Size:	1368 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208668).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138882.1 , NP_620237.1
RefSeq Size:	1743 bp
RefSeq ORF:	1371 bp
Locus ID:	85311
UniProt ID:	P97535
Cytogenetics:	11q21
Gene Summary:	serine specific phospholipase that catalyzes the hydrolysis of a fatty acyl residue at the sn-1 position of lysophosphatidylserine and phosphatidylserine [RGD, Feb 2006]



[View online »](#)