

Product datasheet for **RR208686L4V**

Chid1 (NM_001047854) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Chid1 (NM_001047854) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Chid1
Synonyms:	Cc1-9; RGD1308955
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001047854
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208686).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001047854.2 , NP_001041319.2
RefSeq Size:	1446 bp
RefSeq ORF:	1191 bp
Locus ID:	293628
UniProt ID:	A0JPO9
Cytogenetics:	1q41
Gene Summary:	Saccharide- and LPS-binding protein with possible roles in pathogen sensing and endotoxin neutralization. Ligand-binding specificity relates to the length of the oligosaccharides, with preference for chitotetraose (in vitro) (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)