

## Product datasheet for **RR208717L3V**

### F2r (NM\_012950) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	F2r (NM_012950) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	F2r
Synonyms:	Par1; TRGPC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012950
ORF Size:	1296 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208717).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012950.2</a> , <a href="#">NP_037082.2</a>
RefSeq Size:	3442 bp
RefSeq ORF:	1299 bp
Locus ID:	25439
Cytogenetics:	2q12
Gene Summary:	plays a role in blood coagulation, inflammatory response, and cell proliferation; mediates signal transduction via MAP kinase and other signaling pathways [RGD, Feb 2006]



[View online »](#)