

Product datasheet for RR208717L4

F2r (NM_012950) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	F2r (NM_012950) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	F2r
Synonyms:	Par1; TRGPC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208717).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 25439

Cytogenetics: 2q12

Gene Summary: plays a role in blood coagulation, inflammatory response, and cell proliferation; mediates signal transduction via MAP kinase and other signaling pathways [RGD, Feb 2006]