

Product datasheet for **RR208807L3V**

Mef2d (NM_030860) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mef2d (NM_030860) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Mef2d
Synonyms:	MGC93400
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_030860
ORF Size:	1521 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208807).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030860.2 , NP_110487.2
RefSeq Size:	2165 bp
RefSeq ORF:	1524 bp
Locus ID:	81518
UniProt ID:	O89038
Cytogenetics:	2q34
Gene Summary:	DNA-binding protein; may play a role in regulation of muscle-specific gene expression and differentiation [RGD, Feb 2006]



[View online »](#)