

Product datasheet for **RR208819**

Neil1 (NM_001025754) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Neil1 (NM_001025754) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Neil1
Synonyms: MGC116267
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR208819 representing NM_001025754
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCGGAGGGCCAGAGCTGCACCTGGCCAGCCACTTTGTGAATGAGATATGTAAGGGCTGGTATTCCG
GTGGGTGTGTGGAGAAGTCTCTGTGAGCCGAAACCCGGAGGTGCCCTTTGAGAGCAGTGCCTACCATC
CTCAGCTTTAGCCCGAGGCAAGGAGCTGCGCTTGACATTGAGCCCCTGCCTGGGTCCCATCCCCCTCAG
AAACCACTGTCCCTTGTCTTCCGCTTTGGTATGTCTGGATCCTTCCAGCTGGTACCCGCAGAAGCACTGC
CGCGCCACGCCATCTACGTTTCTATACGGCCCCACCTGCCCCCGGCTTGCCCTTTGCTTCGTAGACAT
CCGCCGCTTTGGCCACTGGGACCCTGGAGGTGAATGGCAGCCAGGCCGCGGACCCTGTGTCTTGTGGAG
TATGAACGGTTCAGAGAGAACGTGCTTCGGAACCTATCGGACAAAGCCTTTGACCGGCCATCTGCGAGG
CCTTGTGGACCAGAGTTCTTCAATGGAATCGGCAACTATCTGCGGGCAGAGATCTTGTATCGGCTGAA
GATCCCCCTTTTGAGAAGGCCGTACAGTTCTAGAGGCCCTGCAACAGTGCCGGCTGAGCCAGAGCTG
ACCCTGAGCCAGAAGGTCAAGGCCAAACTGCAGAACCCAGACCTGCTGGAAGTGTGCTGCTGGTCCCA
AGGAAGTGGTTCAGCTGGGGGGCAAAGGCTATGGCCAGAGCGTGGAGAGGAAGATTTGTGCTGCCCTTCG
AGCCTGGCTTCGGTCTATGCTGTGCCAGGCATGAGCTCCCTGAGAGACCGGCATGGCCGTACCATTTGG
TTCCAGGGTGATCCTGGACCCTTGGCACCCAAAGGGGGCAGATCCCGAAAAAAGAAGTCACAGGAAACAC
AGCTGGAGGCTGAGGACAGAAAAGAGGACCCTCCAAGCAAGTCCATTTCCAGGACTCAGAGGGCCAGGAA
GCACTGTCTAAGAATGTACCTCAACAGTCTGAAGGGGCCAGCCTCCACAAAACCTGGAACCCCTGCA
GTCCCTGAGAAAGGAAAGAGAAGGGGGAGACGGGCGAGCACAGGCCGCCACAGAGGACCAAGAGTAAGC
CTGATACTCAACCAGGGAGGCTGGGGAGATTTACGCCTCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR208819 representing NM_001025754
Red=Cloning site Green=Tags(s)


MPEGPELHLASHFVNEICKGLVFGGCVKSSVSRNPEVPFESSAYHISALARGKELRLTSLPLPGSHPPO
 KPLSLVFRFGMSGFQLVPAEALPRHAHLRFYTAPPAPRLALCFVDIRRFGHWDPPGGEWQPRGRPCVLL
 YERFRENVLRLNSDKAFDRPICEALLDQRFNFNGIGNYLRAEILYRLKIPPFKARTVLEALQQCRLSPEL
 TLSQKVKAKLQNPDLLELCHLVPKEVVQLGGKGYGPERGEEDFAAFRAWLRICYAVPGMSSLRDRHGRTI
 W FQGDPPGLAPKGGRSRKKKSQETQLEADRKEDPPSKSISRTQARKHCPKNVPQQSEGASLPQNLETPA
 VPEKGRRRRASTGRHRGPKSKPDTQPREAGEISAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001025754

ORF Size:	1161 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001025754.1 , NP_001020925.1
RefSeq Size:	1768 bp
RefSeq ORF:	1164 bp
Locus ID:	367090
Cytogenetics:	8q24
MW:	43.3 kDa