

## Product datasheet for **RR208839L4V**

### **Tm6sf2 (NM\_001127654) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Tm6sf2 (NM_001127654) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tm6sf2
Synonyms:	MGC188225
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001127654
ORF Size:	1134 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208839).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001127654.1</a> , <a href="#">NP_001121126.1</a>
RefSeq Size:	1412 bp
RefSeq ORF:	1137 bp
Locus ID:	689029
UniProt ID:	<a href="#">B0BNG2</a>
Cytogenetics:	16p14
Gene Summary:	Regulator of liver fat metabolism influencing triglyceride secretion and hepatic lipid droplet content. May function as sterol isomerase.[UniProtKB/Swiss-Prot Function]



[View online »](#)