

Product datasheet for **RR208950**

Man1a2 (NM_001106452) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Man1a2 (NM_001106452) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Man1a2
Synonyms:	Man1b
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR208950 representing NM_001106452
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACTACCCAGCGCTGCTGCCCTGTCTGGTCGTAGGATACCACCTCTGAACCTGGGGCCGCCCTCT
TCCCACATCACAGGGCTACCTTGAGACTCTCCGAGAAGTTTATCCTTCTCCTCATCCTTAGTGCCTTCAT
CACCTGTGTTTTGGGGCATTCTTCTTCCAGACTCTTCAAAACACAAACGCTTTGATCTGGGCTTA
GAAGATGTGTTAATTCCTCACGTAGATGCGGGCAAAGGGGCTAAAGACCCCGAGTCTTCTGATTCATG
GACCCGACGAACACAGACACAGGAAGAAGAAGAGCGTCTGAGAAATAAGATTAGAGCTGATCATGAGAA
AGCCCTGGAAGAAGCAAAAGAAAAATTAAGAAAGTCAAGAGAGGAAATCCGTGCAGAAATTCAGACAGAG
AAAAATAAAGTAGCCCAAGCAATGAAGACAAAAGAGACCAGGGTACTGCCACCTGTCCTGTCCACAGC
GTGTAGGGGTGAGTGGTGGGACCCAGAAGACGCGGAGATCAGGAAGAAAAGAGACAAAATTAAGAGAT
GATGAAACATGCCTGGGATAATTACAGAACATACGGATGGGGCCATAATGAGCTCAGACCTATTGCAAGG
AAAGGGCATTCCACTAACATATTCGGAAGCTCACAGATGGGGCCACCATAGTGGATGCTTTGGATACCC
TTTATATCATGGGGCTTCATGATGAATTCATGGATGGGCAAAGATGGATTGAAGACAACCTTGATTCAG
TGTGAACTCGGAAGTGTCTGTGTTGAAGTTAACATTCGCTTTATTGGAGGGCTCCTTGCAGCATATTAC
CTGTCAGGAGAGGAAATTTCAAGACTAAAGCAGTGCAGTTGGCTGAGAAGCTCCTTCTGCCTTTAAACA
CTCCTACTGGGATTCCTTGGGCAATGGTGAACCTGAAAAGTGGAGTAGGTCGGAAGTGGGGCTGGGCATC
TGCAGGCAGCAGCATCCTGGCTGAATTTGGCACCTGCACATGGAGTTCGTGCACCTCAGCTACTTGACT
GGTGACCTGACTTACTATAATAAGGTGATGCACATCCGGAACTACTACAGAAAATGGAACGCCAAATG
GTCTTTATCCAAATTTTAAACCCGAGGACAGGGCGCTGGGGTCAGTATCACACATCAGTTGGTGGCCT
AGGTGACAGTTTTTATGAATACTTACTAAAAGCATGGCTGATGTCAGATAAAACAGACCAGGAGCAAGA
AGGATGTACGACGACGCTGTTGAGGCTATAGAGAAACATCTTATTAAGAAGTCCCGAGGAGTCTAGTCT
TTATTGGAGAATGAAAAATGGACACTTGGAAAGGAAGATGGGTCATTTGGCCTGCTTTGCTGGGGAAAT
GTTTGCCTTGGAGCGGATGGTTCCAGAAAGGATAAAGCTGGCCACTATTTAGAAGTGGGGCAGAAAT
GCACGAACATGCCATGAGTCATATGACAGAAGTGCATTAAGTGGTCTGAGTCAATCAAGTTTGATG
GTGCAAGTGAAGCAGTGGCTGTGCGCAAGCTGAGAAGTATTACATCCTTCGTCCGGAAGTATCGAAAC
CTATTGGTATCTATGGCATTACCCACGACCAAGATACAGGCAGTGGGGCTGGGAAGCCGCATTGGCT
ATTGAGAAGTCTGCCGAGTCAGTGGTGGTTTTCTGGTGTCAAGGATGTCTACTCCCGACTCCTGCGC
ACGACGACGTGCAGCAGAGCTTTTTTCTTCCGAGACGTTAAAATACTTGTACCTGCTGTTCTCTGGCGA
TGACCTCTACCTTTAGACCACTGGGTGTTAATACAGAGGCTCATCCTCTGCCGGTGTGCGTTTAGCC
AACACCACTCTTTCAGGTAACCCTGCTGTTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR208950 representing NM_001106452
 Red=Cloning site Green=Tags(s)


MTTPALLPLSGRRIPPLNLGPPSPFHHRATLRLSEKFILLILSAFITLFCGAFFFLPDSSKHKRFDLGL
 EDVLIPHVDAGKGAKDPGVFLIHGPDEHRHREERLRNKRADHEKALEEAKEKLRKSREEIRAEIQTE
 KNKVAQAMKTKETRVLPVPVQVRVGVSGGDPEDAEIRKKRDKIKEMMKHAWDNRYTYGWHNELRPIAR
 KGHSTNIFGSSQMGATIVDALDTLYIMGLHDEFMDGQRWIEDNLDVSVNSEVSVFEVNIREFIGLLAAYY
 LSGEEIFKTKAVQLAEKLLPAFNPTGIPWAMVNLKSGVGRNWGWASAGSILAEFGTLHMEFVHLSYLT
 GDLYTYNKMVHIRKLLQKMERPNGLYPNYLNPRTGRWQYHTSVGGLGDSFYEYLLKAWLMSDKTDHEAR
 RMYDDAVEAIEKHLIKKSRGGLVFIGEWKNGHLERKMGLACFAGGMFALGADGSRKDKAGHYLELGAEI
 ARTCHESYDRALKLGPESFKFDGAVEAVAVRQAEKYYILRPEVIETYWYLWRFTHDPRYRQGWEEAALA
 IEKSCRVSGGFSGVKDVYSPTPAHDDVQQSFFLAETLKYLYLLFSGDDLPLDHWVFNTEAHLPLVRLA
 NTTLSGNPAVR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106452

ORF Size: 1923 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106452.1](#), [NP_001099922.1](#)

RefSeq Size: 3482 bp
RefSeq ORF: 1926 bp
Locus ID: 295319
Cytogenetics: 2q34
MW: 72.8 kDa