

Product datasheet for **RR208985L4V**

Tspan2 (NM_022589) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tspan2 (NM_022589) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tspan2
Synonyms:	Tspan-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022589
ORF Size:	663 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208985).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022589.1 , NP_072111.1
RefSeq Size:	3854 bp
RefSeq ORF:	666 bp
Locus ID:	64521
UniProt ID:	Q9J1W1
Cytogenetics:	2q34
Gene Summary:	may mediate signalling in oligodendrocytes in the early stages of their terminal differentiation into myelin-forming glia; may function in stabilizing the mature sheath [RGD, Feb 2006]



[View online »](#)