

Product datasheet for **RR209053L4V**

Senp7 (NM_001105888) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Senp7 (NM_001105888) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Senp7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001105888
ORF Size:	3111 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209053).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001105888.1 , NP_001099358.1
RefSeq Size:	4876 bp
RefSeq ORF:	3114 bp
Locus ID:	288167
UniProt ID:	D3ZF42
Cytogenetics:	11q12
Gene Summary:	Protease that deconjugates SUMO2 and SUMO3 from targeted proteins, but not SUMO1. Catalyzes the deconjugation of poly-SUMO2 and poly-SUMO3 chains. Has very low efficiency in processing full-length SUMO proteins to their mature forms (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)