

Product datasheet for **RR209070L3V**

Exoc3 (NM_001024964) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Exoc3 (NM_001024964) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Exoc3
Synonyms:	rSec6; Sec6; Sec6l1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001024964
ORF Size:	2265 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209070).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024964.1 , NP_001020135.1
RefSeq Size:	2420 bp
RefSeq ORF:	2268 bp
Locus ID:	252881
UniProt ID:	Q62825
Cytogenetics:	1p11
Gene Summary:	may be involved in the later stages of vesicle trafficking [RGD, Feb 2006]



[View online »](#)