

Product datasheet for **RR209091L4V**

Vdr (NM_017058) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Vdr (NM_017058) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Vdr
Synonyms:	Nr1i1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017058
ORF Size:	1269 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209091).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017058.1 , NP_058754.1
RefSeq Size:	2043 bp
RefSeq ORF:	1272 bp
Locus ID:	24873
UniProt ID:	P13053
Cytogenetics:	7q36
Gene Summary:	binds vitamin D3; acts as a transcriptional regulator of the osteocalcin gene [RGD, Feb 2006]



[View online »](#)