

Product datasheet for **RR209096L3V**

Cep55 (NM_001025646) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cep55 (NM_001025646) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cep55
Synonyms:	RGD1305340
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001025646
ORF Size:	1386 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209096).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025646.1 , NP_001020817.1
RefSeq Size:	2481 bp
RefSeq ORF:	1389 bp
Locus ID:	294074
UniProt ID:	Q4V7C8
Cytogenetics:	1q53



[View online »](#)

Gene Summary:

Plays a role in mitotic exit and cytokinesis. Recruits PDCD6IP and TSG101 to midbody during cytokinesis. Required for successful completion of cytokinesis. Not required for microtubule nucleation. Plays a role in the development of the brain and kidney.[UniProtKB/Swiss-Prot Function]