

## Product datasheet for **RR209128L4V**

### Adh1 (NM\_019286) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Adh1 (NM_019286) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adh1
Synonyms:	Adh; Adh1a; Adh1c
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019286
ORF Size:	1128 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209128).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019286.4</a> , <a href="#">NP_062159.3</a>
RefSeq Size:	1806 bp
RefSeq ORF:	1131 bp
Locus ID:	24172
UniProt ID:	<a href="#">P06757</a>
Cytogenetics:	2q44
Gene Summary:	alpha subunit of class I alcohol dehydrogenase; metabolizes a wide variety of substrates including ethanol, hydroxysteroids and lipid peroxidation products [RGD, Feb 2006]



[View online »](#)