

Product datasheet for **RR209140L4V**

Kcnip4 (NM_181365) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kcnip4 (NM_181365) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kcnip4
Synonyms:	Kchip4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_181365
ORF Size:	750 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209140).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181365.2 , NP_852030.1
RefSeq Size:	1289 bp
RefSeq ORF:	753 bp
Locus ID:	259243
UniProt ID:	Q99MG9
Cytogenetics:	14q11
Gene Summary:	a potassium channel interacting protein which modulates the opening and closing, and disrupts rapid inactivation of the Kv4 potassium channel [RGD, Feb 2006]



[View online »](#)