

## Product datasheet for **RR209142L3V**

### Siah2 (NM\_134457) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Siah2 (NM_134457) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Siah2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134457
ORF Size:	975 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209142).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_134457.1</a> , <a href="#">NP_604452.1</a>
RefSeq Size:	2048 bp
RefSeq ORF:	978 bp
Locus ID:	140593
UniProt ID:	<a href="#">Q8R4T2</a>
Cytogenetics:	2q26
Gene Summary:	a synaptophysin binding protein and also interact with E2 ubiquitin-conjugating enzyme UbcH8; involved in the regulating the ubiquitination and degradation of synaptophysin [RGD, Feb 2006]



[View online »](#)