

Product datasheet for **RR209171L3V**

Nckap1 (NM_031618) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nckap1 (NM_031618) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nckap1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031618
ORF Size:	3384 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209171).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031618.1 , NP_113806.1
RefSeq Size:	4270 bp
RefSeq ORF:	3387 bp
Locus ID:	58823
UniProt ID:	P55161
Cytogenetics:	3q24
Gene Summary:	interacts with Nck and may be involved in the downstream signaling of Nck [RGD, Feb 2006]



[View online »](#)