

## Product datasheet for **RR209173L3V**

### **Bgn (NM\_017087) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bgn (NM_017087) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bgn
Synonyms:	BSPG1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017087
ORF Size:	1107 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209173).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017087.1</a> , <a href="#">NP_058783.1</a>
RefSeq Size:	2446 bp
RefSeq ORF:	1110 bp
Locus ID:	25181
UniProt ID:	<a href="#">P47853</a>
Cytogenetics:	Xq37
Gene Summary:	proteoglycan; may play a role in vascular remodeling [RGD, Feb 2006]



[View online »](#)