

Product datasheet for **RR209202L4V**

Etfb (NM_001004220) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Etfb (NM_001004220) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Etfb
Synonyms:	MGC94808
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001004220
ORF Size:	765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209202).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004220.1 , NP_001004220.1
RefSeq Size:	867 bp
RefSeq ORF:	768 bp
Locus ID:	292845
UniProt ID:	Q68FU3
Cytogenetics:	1q22
Gene Summary:	human homolog plays a role in fatty acid metabolism and beta-oxidation [RGD, Feb 2006]



[View online »](#)