

## Product datasheet for **RR209223**

### Gltp (NM\_001134413) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gltp (NM\_001134413) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Gltp  
**Synonyms:** RGD1564442  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR209223 representing NM\_001134413  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGCGCTGCTGGCCGAGCACCTGCTCAAGCCGCTACCCGAGACCGGCAGATCGAGACCGGCCCTTCT  
TGGAGGCGGTGGCCACCTGCCGCCCTTCTCGATTGCCTGGGTCCCCGGTGTTCACGCCCATCAAGGC  
AGACATAAGCGGTAACATCACAAAATCAAAGCGTCTATGACTGACCCGGCCAAGTCAAGACTCTG  
CAGAACATCCTGGAGGTGAGAAGGGGATGTATGGGGCGGAGTGCCCAAAGTGGGGGCCACGCTGGCAC  
TGCTGTGGCTGAAAAGGGGCTGCGCTTCCAGGTGTTCTGCAGAGCATCTGCGATGGCGAACGGGA  
CGAGAACCACCCCAACCTCATCCGGTAAACGCCAACAAGGCCTATGAGATGGCGCTAAAGAAGTATCAC  
GGCTGGCTGGTGAGAAGATCTCAAGGCCGACTGTATGCAGCACCATAACAAGTCTGACTTCTAAAGG  
CGCTCTCAAGGGGCAGAACGTGACGGAGGAGGAGTGCTGGAGAAGATCCGCCTTCTCTGGTCAACTA  
CACGGCCACCATCGACGCCATCTATGAAATGTATACCAAGATGAACGCGGAGCTTGACTACACCGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR209223 representing NM\_001134413  
Red=Cloning site Green=Tags(s)

MALLAEHLLKPLPADRQIETGPFLEAVAHLPFFDCLGSPVFTPIKADISGNITKIKAVYDTPAKFKTL  
QNILEVEKMGYAEWPKVGATLALLWLRGLRFIQVFLQSIDGERDENHPNLIRVNANKAYEMALKKYH  
GWLVQKIFKAALYAAPYKSDFLKALSKGQNVTEEECLEKIRLFLVNYTATIDAIYEMYTKMNAELDYTV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI




[View online »](#)

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001134413

**ORF Size:** 627 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001134413.1](#), [NP\\_001127885.1](#)

**RefSeq Size:** 1666 bp

**RefSeq ORF:** 630 bp

**Locus ID:** 288707

**UniProt ID:** [B0BNM9](#)

**Cytogenetics:** 12q16

**MW:** 23.7 kDa

**Gene Summary:** Accelerates the intermembrane transfer of various glycolipids. Catalyzes the transfer of various glycosphingolipids between membranes but does not catalyze the transfer of phospholipids. May be involved in the intracellular translocation of glucosylceramides (By similarity).[UniProtKB/Swiss-Prot Function]