

## Product datasheet for **RR209232**

### Rpf2 (NM\_001106391) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rpf2 (NM\_001106391) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Rpf2  
**Synonyms:** Bxdc1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR209232 representing NM\_001106391  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGATGCTTTGGACAAAGTCCTAAAGCCTAAAACAAAAAGGGCTAAGAGATTCCTTGAGAAGAGAGAAC  
CAAACTCAGTGAGAACATTAAGAATGCCATGCTGATTAAGGCGGGAATGCGAATGCAACGGTGACGCA  
AGTACTGAGAGACGTGTATGCACTGAAGAAACCGTATGGAGTTCTCTATAAAAAGAAAAACATTACAAGA  
CCATTTGAGGATCAGACGTCCTGGAGTTCTTTCAAAGAAGTCGGACTGTTCTCTGTTTCATGTTTGGCT  
CCCACAATAAGAAGAGGCCCAATAATCTCGTGATAGGTCGTATGTATGACTACCATGTGCTGGATATGAT  
TGAAGTAGGGATTGAGAAATTCGTGCTCTGAAAGATATTAAGACAAGTAAGTGCCTGAGGGAACAAAAG  
CCCATGTTGATATTTGCTGGTGATGATTTTGATGTAACAGAAGACTACAGAAGACTGAAAAATCTTCTTA  
TTGATTTCTTCAGAGGCCCCACAGTGTCAAACGTCCTGGCTGGATTAGAGTATGTTCTGCACTTAC  
GGCACTGAATGGGAAGGTTTACTTTCAAGCTATAAGCTGCTGCTGAAGAAATCTGGGTGCAGAACCCG  
AGAATCGAGCTAGAAGAGATGGGCCATCGTTGGATCTGGTTGTGAGGAGGACACACTTAGCATCGGACG  
ACCTTTATAAAGTGTCCATGAAAGTGCCCAAAGCTCTCAAGCCGAAGAAGAGGAAGAATATCTCTCAGGA  
TACTTTTGGTACAACGTTTGAAGGATTCATATGCAGAAGCAAGACCTAAGCAAACACTACAGACGAGGAAA  
GTGAAGGGTTGAAGAAGCGCCTGCAGAGAATGGAGTAGAAGACCAGGGGAAGAATCAAAGAGAAATTA  
AGAAAAGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR209232 representing NM\_001106391  
 Red=Cloning site Green=Tags(s)

MDALDKVLKPKTKRAKRFLEKREPKLSENIKNAMLIKGGNANATVTQVLRDVYALKKPYGVLYKKKNITR  
 PFEDQTSLEFFSKKSDCSLFMFGSHNKKRPNNLVIGRMYDYHVLDMIELGIEKFVSLKDIKTSKCPEGTK  
 PMLIFAGDDFDVTEYRRLKNLLIDFFRGPTVSNVRLAGLEYVLHFTALNGKVYFRSYKLLKSGCRTP  
 RIELEEMGPSLDLVVRRTHLASDDLYKLSMKVPKALKPKKRKNI SQDTFGTTFGRIHMQKQDL SKLQTRK  
 VKGLKKRPAENGVEDQGKKSRIKKS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001106391

ORF Size: 918 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001106391.2</a> , <a href="#">NP_001099861.1</a>
<b>RefSeq Size:</b>	1226 bp
<b>RefSeq ORF:</b>	921 bp
<b>Locus ID:</b>	294436
<b>Cytogenetics:</b>	20q12
<b>MW:</b>	35.3 kDa