

Product datasheet for **RR209245L4V**

Dnajib9 (NM_012699) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dnajib9 (NM_012699) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dnajib9
Synonyms:	ERdj4; Mdg1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012699
ORF Size:	666 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209245).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012699.2 , NP_036831.2
RefSeq Size:	1936 bp
RefSeq ORF:	669 bp
Locus ID:	24908
UniProt ID:	P97554
Cytogenetics:	6q21
Gene Summary:	mammalian chaperone; involved in the control of cell cycle arrest taking place during terminal cell differentiation and under stress conditions [RGD, Feb 2006]



[View online »](#)