

Product datasheet for **RR209257**

Borcs5 (NM_001108650) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Borcs5 (NM_001108650) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Borcs5
Synonyms: Loh12cr1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR209257 representing NM_001108650
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGCAGCGAGCAGAGCGGGAGGCCGAGAGCCGCCCGGCGATCTGAACGCCTCAGTGACTCCATCTC
CGGCTAAACATCGGGCAAGATGGATGATATCGTGGTGGTAGCTCAGGGCTCTCAGGCCTCAAGGAATG
CAGCAATGATCCCGATGTCATCAAGCTGCAAGAGATTCCAACCTCCAGCCCCTTTAAAAGGGCTCTTG
AGTGGACAGAGCTCCCAACAAACGCAAGGCTGGAGAAGCTGGACTCTCAGCAGGTGCTGCAGCTCTGCC
TCCGGTACCAAGACCACCTGCACCAGTGCAGAGGCCGTCGCTTTCGACCAGAATGCTCTGGTCAAGCG
AATCAAAGAGATGGATCTGTCCGTGAAACCTCTTTGCTTCATGCAGGAGCGCCAGAAAAGATACGCC
AAATACGCTGAACAGATCCAGAAGGTGAACGAGATGTCGGCCATCCTGCGCCGTATCCAGATGGGCATAG
ACCAGACGGTGCCACTCATGGAGAGGCTCAACAGCATGCTGCCAGAGGCCGAGCGCCTGGAGCCCTTCAG
CATGAGGCCTGAGCGGGAGCGCCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR209257 representing NM_001108650
Red=Cloning site Green=Tags(s)

MGSEQSAEASRPGLNASVTPSPAKHRAKMDDIVVVAQGSQASRNVSNPDVVIKLEIPTFQPLLKGLL
SGQTSPNARLEKLDSSQVQLQLCLRYQDHLHQCAEAVAFDQNALVKRIKEMDLSVETLFCFMQERQKRYA
KYAEQIQKVNEMSAILRRIQMGIDQTVPLMERLNSMLPEAERLEPFMRPERERH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108650

ORF Size: 585 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001108650.2](#), [NP_001102120.1](#)

RefSeq Size: 1852 bp

RefSeq ORF: 588 bp

Locus ID: 362452

Cytogenetics: 4q43

MW: 22.2 kDa