

Product datasheet for **RR209262**

TPD52I2 (NM_198744) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tpd52I2 (NM_198744) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tpd52I2
Synonyms: MGC73020
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR209262 representing NM_198744
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATTCTGCTAGCCAAGATATCAACCTGAATTCCTAACAAGGTGTGCTGTCTGATTTTATGACTG
ATGTCCTGTTGACCCAGGTGTGGTCCACCGGACTCCAGCTGTAGAGGGCCTGACCGAGGTGGAGGAAGA
AGAGCTCCGGGCTGAGCTTGCTAAGGTGGAAGAGGAAATTGTCACTCTGCGCCAGGTGCTGGCAGCCAAA
GAGAGGCACTGTGGAGAGCTGAAACGGAGACTGGGCCTCTCCACATTAGGGGAGCTGAAGCAGAACCTGT
CTAGGAGCTGGCATGATGTGCAGGGCTCTACTGCCTATGTGAAAACGCTCTGAGAACTTGAGAGTGGAA
TGAGAAAGTGACGCAGTCTGACCTCTACAAGAAGACTCAGGAACTCTTTCACAGGCAGGGCAGAAAACA
TCAGCTGCTCTGTCCACCATGGGCTCTGCCATCAGCAGGAAGCTCGGAGATATGAGCAGCTACTCCATCC
GACACTCGATAAGTATGCCTGTCATGAGGAACTCTGCAACCTTCAAGTCATTTGAAGATCGAGTGGGGAC
CATAAAGTCTAAGGTTGTGGCGGCAGAGAGAATGGCAGCGATACCCTCCCTTCTCACCTGGAAGTGGT
GACCAGACATTGCCGGATCATGCGCCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR209262 representing NM_198744
Red=Cloning site Green=Tags(s)

MDSASQDINLNSPNKGVLSDFMTDVPVDPGVVHRTPAVEGLTEVEEEELRAELAKVEEEIVTLRQVLAAK
 ERHCGELKRRLGLSTLGELKQNLRSWHDVQGSTAYVKTSEKLGEWNEKVTQSDLYKKTQETLSQAGQKT
 SAALSTMGSAISRKLGDMSYSIRHSISMPVMRNSATFKSFEDRVGTIKSKVVGRENGSDTLPPSSPGSG
 DQTLPDHAPF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_198744

ORF Size: 660 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198744.2 , NP_942039.1
RefSeq Size:	1407 bp
RefSeq ORF:	663 bp
Locus ID:	296480
UniProt ID:	Q6PCT3
Cytogenetics:	3q43
MW:	24 kDa
Gene Summary:	human homolog is a member of a family of proteins that may play a role in cell proliferation [RGD, Feb 2006]