

Product datasheet for **RR209281L4V**

Napepld (NM_199381) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Napepld (NM_199381) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Napepld
Synonyms:	NAPE-PLD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199381
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209281).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199381.1
RefSeq Size:	1191 bp
RefSeq ORF:	1191 bp
Locus ID:	296757
UniProt ID:	Q769K2
Cytogenetics:	4q11
Gene Summary:	an atypical phospholipase D that catalyzes the conversion of N-acyl-phosphatidylethanolamines to anandamide and other N-acylethanolamines but has no activity with phosphatidylcholine or phosphatidylethanolamine [RGD, Feb 2006]



[View online »](#)