

## Product datasheet for **RR209313L4V**

### Gjb1 (NM\_017251) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gjb1 (NM_017251) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gjb1
Synonyms:	Cx32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017251
ORF Size:	849 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209313).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017251.2</a> , <a href="#">NP_058947.1</a>
RefSeq Size:	1555 bp
RefSeq ORF:	852 bp
Locus ID:	29584
UniProt ID:	<a href="#">P08033</a>
Cytogenetics:	Xq22
Gene Summary:	component of gap junctions, which connect the cytoplasm of adjacent cells; involved in intercellular transport of small hydrophilic molecules; facilitates cell-cell communication [RGD, Feb 2006]



[View online »](#)