

Product datasheet for **RR209344L3V**

Rcbtb2 (NM_199084) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rcbtb2 (NM_199084) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rcbtb2
Synonyms:	Chc1l
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_199084
ORF Size:	1653 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209344).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199084.1 , NP_954515.1
RefSeq Size:	1971 bp
RefSeq ORF:	1656 bp
Locus ID:	290363
UniProt ID:	Q6P798
Cytogenetics:	15q11
Gene Summary:	human homolog may be a guanine nucleotide exchange factor related to the regulator of chromosome condensation (RCC1) gene [RGD, Feb 2006]



[View online »](#)