

## Product datasheet for **RR209389**

### Ube2k (NM\_001106006) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ube2k (NM\_001106006) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ube2k  
**Synonyms:** E2-25k; Hip2; Hypg; Lig  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR209389 representing NM\_001106006  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCCAACATCGCGGTGCAGCGAATCAAGCGGGAGTTCAAGGAGGTGCTGAAGAGCGAGGAGACGAGCA  
 AAAATCAAATTAAGTAGATCTTGTAGATGAGAATTTACAGAATTAAGAGGAGAAATAGCAGGACTCC  
 AGACACACCGTATGAAGGTGGAAGATATCAACTAGAGATAAAAAACCAGAAACATATCCATTTAACCC  
 CCTAAGGTCCGGTTTATCACTAAAATATGGCACCTAATATTAGTCCGTACAGGGGCTATTTGTTGG  
 ATATCCTGAAAGATCAATGGGCAGCAGCAATGACTCTGCGCACAGTATTATTGTCATTGCAAGCACTGCT  
 GGCCGCTGCAGAACCAGATGACCCCAAGATGCAGTAGTACAAATCAGTACAAACAAAATCCTGAAATG  
 TTCAAGCAGACAGCTCGACTTTGGGCACACGTGTATGCTGGAGCACCAGTGTCTAGTCCAGAATACACCA  
 AAAAAATAGAAAACCTCTGTGCTATGGGTTTTGATAGGAACGCAGTAATAGTGGCCTTGTCTTCAAATC  
 ATGGGATGTAGAGACTGCAACAGAACTGCTTCTGAGTAAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR209389 representing NM\_001106006  
 Red=Cloning site Green=Tags(s)

MANIAVQRIKREFKEVLKSEETSKNQIKVDLVDFTELRLGEIAGPPDTPYEGGRYQLEIKIPETYPFNP  
 PKVRFITKIWHPNISSVTGAICLDILKDQWAAAMTLRTVLLSLQALLAAEPDDPQDAVVANQYKQNP  
 FKQTARLWAHVYAGAPVSSPEYTKKIENLCAMGFDRNAVIVALSSKSWDVETATELLLSN


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_001106006

**ORF Size:** 600 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001106006.1](#), [NP\\_001099476.1](#)

**RefSeq Size:** 4971 bp

**RefSeq ORF:** 603 bp

**Locus ID:** 289623

**Cytogenetics:** 14p11

**MW:** 22.4 kDa