

Product datasheet for **RR209394L3V**

Ccdc80 (NM_022543) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccdc80 (NM_022543) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ccdc80
Synonyms:	CL2; Dro1; Ssg1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022543
ORF Size:	2847 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209394).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022543.2 , NP_071988.2
RefSeq Size:	4114 bp
RefSeq ORF:	2850 bp
Locus ID:	64387
UniProt ID:	Q6QD51
Cytogenetics:	11q21
Gene Summary:	estrogen regulated gene; overexpression correlates with mammary carcinogenesis [RGD, Feb 2006]



[View online »](#)