

## Product datasheet for **RR209413L3V**

### Ada (NM\_130399) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ada (NM_130399) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ada
Synonyms:	MGC108610
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_130399
ORF Size:	1056 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209413).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_130399.2</a> , <a href="#">NP_569083.1</a>
RefSeq Size:	1417 bp
RefSeq ORF:	1059 bp
Locus ID:	24165
UniProt ID:	<a href="#">Q920P6</a>
Cytogenetics:	3q42
Gene Summary:	catalyzes the conversion of adenosine to inosine in adenosine metabolism; may act as a neuromodulator in sleep-wake regulation [RGD, Feb 2006]



[View online »](#)