

Product datasheet for **RR209455L4V**

Mrgprf (NM_153722) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Mrgprf
Synonyms:	MrgF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153722
ORF Size:	1029 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RR209455).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_153722.2 , NP_714944.1
RefSeq Size:	2071 bp
RefSeq ORF:	1032 bp
Locus ID:	266762
UniProt ID:	P23749
Cytogenetics:	1q42



Gene Summary: seven transmembrane domain G-protein coupled receptor; most often found in sensory neurons [RGD, Feb 2006]