

Product datasheet for **RR209456L3V**

Gpn3 (NM_201991) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gpn3 (NM_201991) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gpn3
Synonyms:	Atpbd1c
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_201991
ORF Size:	852 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209456).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_201991.1 , NP_973720.1
RefSeq Size:	1176 bp
RefSeq ORF:	855 bp
Locus ID:	360810
UniProt ID:	Q6R518
Cytogenetics:	12q16
Gene Summary:	Small GTPase required for proper localization of RNA polymerase II (RNAPII). May act at an RNAP assembly step prior to nuclear import.[UniProtKB/Swiss-Prot Function]


[View online »](#)