

Product datasheet for **RR209460L4V**

Vezt (NM_001006984) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Vezt (NM_001006984) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Vezt
Synonyms:	MGC94969
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001006984
ORF Size:	2346 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209460).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001006984.1 , NP_001006985.1
RefSeq Size:	2906 bp
RefSeq ORF:	2349 bp
Locus ID:	299738
UniProt ID:	Q5XI52
Cytogenetics:	7q13
Gene Summary:	Plays a pivotal role in the establishment of adherens junctions and their maintenance in adult life. Required for morphogenesis of the preimplantation embryo, and for the implantation process.[UniProtKB/Swiss-Prot Function]


[View online »](#)