

## Product datasheet for **RR209471**

### **Ppil1 (NM\_001034188) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ppil1 (NM\_001034188) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ppil1  
**Synonyms:** MGC116359  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR209471 representing NM\_001034188  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGCATCCCCCAGACACTTGGCAGCCGCCAACGTCTACCTGGAGACTAGCATGGGAATCATTG  
TGCTGGAGCTGTACTGGAAGCACGCTCCAAAGACCTGCAAGAACTTCGCGGAGCTGGCTCGGAGGGGCTA  
CTACAACGGCACCAAGTCCACAGGATCATCAAGGACTTCATGATCCAAGCGGAGACCCGACAGGCACA  
GGTCGAGGTGGCGCCTCTATTTATGGCAAACAGTTTGAAGATGAGCTTCATCCAGACCTGAAGTTCACTG  
GGGCTGGGATTCTTGAATGGCCAATGCAGGACCAGACACCAACGGCAGCCAGTTCTTCGTGACCTTGGC  
CCCCACACAATGGCTGGACGGCAAACATACCATTTTTGGCCGCGTGTGCCAGGGCATAGGGATGGTGAAT  
CGAGTGGGCATGGTGAAACAACTCCCAGGACCGCCTGTGGATGACGTGAAGATCCTGAAGGCTTACC  
CTTCTGGG

**ACGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR209471 representing NM\_001034188  
**Red=Cloning site Green=Tags(s)**

MAAIPPDTWQPPNVYLETSMGIIVLELYWKHAPKTKNFAELARRGYNGTKFHRIKDFMIQGGDPTGT  
GRGGASIYGKQFEDELHPDLKFTGAGILAMANAGPDTNGSQFFVTLAPTQWLDGKHTIFGRVCQIGIMVN  
RVGMVETNSQDRPVDDVKILKAYPSG

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:** SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM\_001034188

ORF Size: 498 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001034188.1</a></u> , <u><a href="#">NP_001029360.1</a></u>
<b>RefSeq Size:</b>	1325 bp
<b>RefSeq ORF:</b>	501 bp
<b>Locus ID:</b>	309651
<b>Cytogenetics:</b>	20p12
<b>MW:</b>	18.3 kDa