

Product datasheet for **RR209474L4V**

Cycs (NM_012839) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cycs (NM_012839) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cycs
Synonyms:	CYCSA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012839
ORF Size:	315 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209474).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012839.2 , NP_036971.1
RefSeq Size:	1228 bp
RefSeq ORF:	318 bp
Locus ID:	25309
UniProt ID:	P62898
Cytogenetics:	4q24
Gene Summary:	a component of the electron transport chain in mitochondria, may function in apoptosis [RGD, Feb 2006]


[View online »](#)