

## Product datasheet for **RR209487L4V**

### **Ccr9 (NM\_172329) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ccr9 (NM_172329) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ccr9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172329
ORF Size:	1107 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209487).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172329.1</a> , <a href="#">NP_758832.1</a>
RefSeq Size:	2568 bp
RefSeq ORF:	1110 bp
Locus ID:	282832
Cytogenetics:	8q32
Gene Summary:	mouse homolog binds CCL25 and is involved in T lymphocyte chemotaxis [RGD, Feb 2006]



[View online »](#)