

Product datasheet for **RR209501L4V**

Crmp1 (NM_012932) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Crmp1 (NM_012932) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Crmp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012932
ORF Size:	1716 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209501).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012932.2 , NP_037064.1
RefSeq Size:	2846 bp
RefSeq ORF:	1719 bp
Locus ID:	25415
UniProt ID:	Q62950
Cytogenetics:	14q21
Gene Summary:	protein that may play a role in neuronal plasticity by transduction of signals from semaphorins [RGD, Feb 2006]



[View online »](#)