

## Product datasheet for **RR209505L3V**

### Mat1a (NM\_012860) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mat1a (NM_012860) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Mat1a
Synonyms:	AdoMet; SADE; SAS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012860
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209505).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012860.2</a> , <a href="#">NP_036992.2</a>
RefSeq Size:	3056 bp
RefSeq ORF:	1191 bp
Locus ID:	25331
UniProt ID:	<a href="#">P13444</a>
Cytogenetics:	16p14
Gene Summary:	enzyme that catalyzes the formation of S-adenosylmethionine [RGD, Feb 2006]



[View online »](#)