

## Product datasheet for **RR209519L4V**

### Cops3 (NM\_001004200) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cops3 (NM_001004200) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cops3
Synonyms:	MGC94156
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001004200
ORF Size:	1269 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209519).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001004200.1</a> , <a href="#">NP_001004200.1</a>
RefSeq Size:	1596 bp
RefSeq ORF:	1272 bp
Locus ID:	287367
UniProt ID:	<a href="#">Q68FW9</a>
Cytogenetics:	10q22
Gene Summary:	subunit of the COP9 signalosome complex; mouse homolog plays a role in early embryonic development [RGD, Feb 2006]



[View online »](#)