

Product datasheet for **RR209524**

Nudt6 (NM_181363) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nudt6 (NM_181363) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nudt6
Synonyms: Gfg
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR209524 representing NM_181363
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGGTGGGCGAGCAGAGCGCGATGGTTGTTCTCGGCTTTGCTCGACGTGGGCGGAGTCGGACTCCGGG
CTAGGCGGCGCACTGCCTCGAGCGGCCTGGAGATCACGGGCTCCTGCGGGGGCGAGCTACAAGGCGAGCT
GGACAGATTTGGGGGCATCTCGGTGCACCTGTGCGGGCACCGCACCTGCACCGGCTGGATGCCGCCGCC
TTCCGGAGGCTCTTGCAGGCCGCCATTCAGCAGTGGCGAGCGGAAGGAAGGATCGCTGCATGGCTGCACA
TCCCCATCCTGCAGAGCCACTTCATCGCCCTGCTGCCTCTCTGGGCTTCTGCTTCCACCACGCAGAACC
GCATTGTCAACGCTGACTCTGTGGTGGGAGAAGGACCCAGCAGACTCCCGGGGTATGCCACACACCAA
GTAGGGGTTGCAGGTGCCGATTTGATGTTAGCACAGGAAAGTACTGGTTGTACAAGATCGAAATAAAT
TGAAAAACATGTGGAAGTTTCCAGGAGGCTGTACAGCCTGGAGAAGATATTGGAGACACAGCAGTCCG
AGAGGTGTTTGAAGAGACTGGTGTCAAGTCAGAATTCAGGTCTCTGTTGAGCATCCGGCAGCAACACAGG
AGCCCTGGGGCCTTTGGGATGTCAGACATGTACCTGATCTGCCGCTGCAGCCGCGTTCCTTACCATCA
ACTTCTGCCAGCAGGAATGCTTGAAGTGTGAATGGATGGATCTAGAAAGCCTGGCCAGGACTAAACACAC
AACCCCATCACAGCAGGGTGGCTAGGCTACTACTATACGGACACAGGAAGGGTTTGACAAGATTGAC
CTCAGCATGGAGAACTCCCTGCAGTATACACAGGCTGTTCTACAAGCTTACCACAGGGGACTGCCTG
AGAGGTACAAGGCTGAAATGGGAACAGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR209524 representing NM_181363
 Red=Cloning site Green=Tags(s)

MWASRARWLF SALLDVGGVGLRARRRTASSGLEITGSCGGELQGELDRFGGISVHLSRHRTLHRLDAAA
 FRRLAQAAIQWRAEGRIAAWLHIPILQSHFIAPAASLGFCHHAEPHLSTLTLWLGEGPSRLPGYATHQ
 VGVAGAVFDVSTRKVLVQDRNKLKNMWFPPGGLSEPGEDIGDTAVREVFEETGVKSEFRSLLSIRQQR
 SPGAFGMSDMYLICRLQPRSFTINFCQQECLKCEWMDLESLARTKHTTPTISRVARLLLYGHREGFDKID
 LSMEELPAVYTGLFYKLYHRGLPERYKAEMGTD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_181363

ORF Size: 939 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_181363.2 , NP_852028.1
RefSeq Size:	1099 bp
RefSeq ORF:	942 bp
Locus ID:	207120
UniProt ID:	P70563
Cytogenetics:	2q25
MW:	35.3 kDa
Gene Summary:	plays a role in the post-transcriptional regulation of FGF-2 [RGD, Feb 2006]